

 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

**Purpose:**

The purpose of this evaluation is to determine selected performance characteristics of the TopOx single-use oxygen mask. The disposable mask features an attached bottle/bag fixture that stores a volume of delivered continuous flow oxygen (CFO) during the user’s exhalation phase. Application of the mask is intended to provide greater FiO<sub>2</sub>% than standard CFO by allowing the user to inhale both the stored oxygen and inspiratory CFO. This evaluation compared FiO<sub>2</sub>% in application of the TopOx with and without the oxygen storage bottle attached in various test conditions. In addition, the effect of mask application on ETCO<sub>2</sub>% was measured as well as any inspiratory restrictions placed on the user due to their tidal volume exhausting the storage volume of the bottle/bag.

**Equipment:**

**Units Under Test (UUT):**

TopOx Single-Use Oxygen Mask – TopOx Ltd, Gloucester, UK



 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

### **Additional Equipment:**

Series 1101 Breathing Simulator- Hans Rudolph, Inc.

Oxygen Analyzer, Model 570A- Servomex

BTC-II Miniature Diaphragm Pump- Hargraves Technology Corp.

Model 4140 Mass Flowmeter- TSI, Inc.

Mannequin Head

Modified with 22mm adapter placed in oro-nasal area

Model 7100 CO<sub>2</sub>SMO ETCO<sub>2</sub>/SpO<sub>2</sub> Monitor- Novamatrix Medical Systems

100% O<sub>2</sub>; 100% CO<sub>2</sub>- Airgas Puritan Bennett

Corr-a-Flex 22mm tubing- Hudson RCI

### **Test Procedures:**

Prior to testing, the TopOx mask was inspected for any product deficiencies or damage; no issues were found.

Using Annex F of the ISO 17510-2 (2007): “Sleep apnoea breathing therapy—Part 2: Masks and application accessories” standard as reference for test setup of a simulated patient, a length of 22mm diameter Corr-a-Flex was cut to measure 140mL volume and used to connect the breathing simulator connection port to the 22mm adapter on the mannequin head.


The TopOx mask was placed on the mannequin head with the headgear securing the mask to the mannequin face. No attempt was made to ensure there was no unintentional leak other than adjusting the fit of the mask to best fit the face.



 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

The Servomex oxygen analyzer was calibrated with room air (20.9%) and 100% oxygen. The ETCO<sub>2</sub> monitor was calibrated using the unit's self-calibration feature.

### Setup and Procedures for FiO<sub>2</sub> Test



Silicone tubing was used to connect the breathing simulator sampling port to the pump inlet port, and also to connect the pump outlet port to the Servomex oxygen analyzer. The pump takes gas samples from the lung and directs them to the oxygen analyzer at a rate of approximately 0.725 LPM.

 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

Standard 25' oxygen tubing was used to connect the 100% O<sub>2</sub> gas source (a cylinder fitted with a low flow regulator) to the flowmeter. Silicone tubing was used to connect the flowmeter to the oxygen inlet port of the mask. Oxygen delivery was regulated to 3 LPM continuous flow (+/- 0.050 LPM), and verified by observing the reading on the flowmeter, adjusting the regulator as necessary.

Three patient conditions were separately tested. Breathing simulator settings are shown below:


	<b>Breathing Simulator Settings</b>				
	<b>Test 1</b>		<b>Test 2</b>		<b>Test 3</b>
<b>Resistance</b>	5		5		5
<b>Compliance</b>	80		80		80
<b>Breath Rate</b>	15		30		40
<b>Amplitude</b>	29.5		34		34.5
<b>Effort Slope</b>	3		3		3
<b>% Inhale</b>	34		50		50
	<b>15 BPM</b>		<b>30 BPM</b>		<b>40 BPM</b>
	<b>500 Vt</b>		<b>400 Vt</b>		<b>300 Vt</b>

For a given test and with the simulator actively breathing, the TopOx was applied and 3 LPM CFO was delivered to the mask. Readings from the Servomex oxygen analyzer were observed until the displayed reading stabilized, at which time the value was recorded.

Following the above tests, each scenario was retested without the TopOx storage bag/bottle attached to the mask.

 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

### Setup and Procedures for ETCO<sub>2</sub> Test



Test setup for ETCO<sub>2</sub> was identical to the FiO<sub>2</sub> test, with two exceptions: The oxygen analyzer and pump were replaced by a CO<sub>2</sub> gas source and regulator, and an ETCO<sub>2</sub> monitor was added to the setup. The CO<sub>2</sub> regulator was connected to the simulator gas sample port via silicone tubing. The ETCO<sub>2</sub> monitor probe was placed between the Corra-Flex tubing and the 22mm adapter on the mannequin head (the probe cannot be seen in the image above, as it is inside the hollow mannequin head).

 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

The breathing simulator was set to run a sinusoidal-type wave pattern at 15bpm with a tidal volume of ~500mL and an I:E ratio of 1:2. With the TopOx mask removed from the mannequin head, gas flow from the CO<sub>2</sub> tank was initiated and ET<sub>CO</sub><sub>2</sub> values were monitored on the capnograph. CO<sub>2</sub> gas flow was adjusted until a ~5 % baseline ET<sub>CO</sub><sub>2</sub> was achieved. This setup was allowed to run for five minutes to ensure baseline ET<sub>CO</sub><sub>2</sub> was stable, at which time baseline ET<sub>CO</sub><sub>2</sub> was recorded.

The TopOx mask was then reapplied to the mannequin head, and the oxygen regulator was adjusted to deliver 3LPM CFO. ET<sub>CO</sub><sub>2</sub> values were monitored on the capnograph. Five minutes were allowed to pass to ensure ET<sub>CO</sub><sub>2</sub> was stable, at which time ET<sub>CO</sub><sub>2</sub> was recorded.

### **Additional Observational Information**

Using the test setup described for the FiO<sub>2</sub> test at 15BPM and 30 BPM, the simulator's Amplitude setting was adjusted upward (this increases tidal volume) until it was noted that the oxygen storage bag of the TopOx was fully collapsed due to the vacuum caused during inhalation. The on-screen graphics of the breathing simulator were monitored to note if there was any added resistance due to the collapse of the bag (identifiable by a sudden reduction or ceasing of flow in the patient flow curve and/or an increase in peak-to-peak airway pressures).

Using the test setup described for the FiO<sub>2</sub> tests, substituting a simulated nose fitted with a Hudson 1104 nasal cannula for the TopOx mask/dummy head fixture, oxygen flow was adjusted through the system so that FiO<sub>2</sub> values matched that recorded from each tested condition with the TopOx mask at 3 LPM CFO. The O<sub>2</sub> flow reported on the flowmeter was recorded to note the liter flow equivalency to the TopOx FiO<sub>2</sub> results.


 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

## Results and Observations:

### FiO<sub>2</sub> Tests

TopOx FiO <sub>2</sub> % at 3LPM CFO				
	Test 1	Test 2	Test 3	
<b>w/ storage</b>	56.6	47.8	46.8	<b>w/ storage</b>
<b>w/o storage</b>	32.8	32.4	33.1	<b>w/o storage</b>
	<b>15 BPM</b>	<b>30 BPM</b>	<b>40 BPM</b>	
	<b>500 Vt</b>	<b>400 Vt</b>	<b>300 Vt</b>	

### FiO<sub>2</sub>% with TopOx Applied, 3 LPM 100% O<sub>2</sub> CFO



 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

**ETCO<sub>2</sub> Test**

		<b>CO2 Rebreathing</b>		
		<b>Baseline ETCO<sub>2</sub>%</b>	<b>Mask ETCO<sub>2</sub>%</b>	<b>% Increase</b>
<b>TopOx</b>	<b>3LPM CFO</b>	5.0	6.4	28

Application of the mask in the setup increased ETCO<sub>2</sub>% from 5.0% baseline to 6.4% (a 28% increase), which is above the limit set by the ISO 17510-2 standard for CPAP masks with pressure applied (20% or less), but is much less than the "fault condition" limit of 60%, where no pressure is applied. While there was an increase in ETCO<sub>2</sub>% due in part to a lack of flushing of the deadspace during exhalation, the one way valves around the mask force volume to be exhaled outside the mask and the mask appears to have enough venting capability that CO<sub>2</sub> is not successively retained.

**Liter Flow Equivalency**

Liter flow equivalency (continuous flow O<sub>2</sub>) needed to achieve the FiO<sub>2</sub> values recorded with the TopOx mask at 3 LPM were as follows:

		<b>Liter Flow Equivalency</b>				
		<b>Test 1</b>	<b>Test 2</b>	<b>Test 3</b>		
<b>FiO<sub>2</sub>%</b>	<b>56.6</b>	<b>47.8</b>	<b>46.8</b>	<b>FiO<sub>2</sub>%</b>		
<b>TopOx</b>	3.0 LPM	3.0 LPM	3.0 LPM	<b>TopOx</b>		
<b>CFO</b>	8.9 LPM	7.5 LPM	6.9 LPM	<b>CFO</b>		
		<b>15 BPM</b>	<b>30 BPM</b>	<b>40 BPM</b>		
		<b>500 Vt</b>	<b>400 Vt</b>	<b>300 Vt</b>		

 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

## Discussion/Conclusions:

Application of the TopOx mask/bottle with 3 LPM CFO resulted in significantly greater  $FiO_{2s}$  than  $FiO_{2s}$  seen with standard 3 LPM CFO. These results indicate that the storage bottle attachment does retain a volume of delivered  $O_2$  during exhalation and that some or all of this volume is inhaled on the subsequent breath, supplementing the 3 LPM CFO that is delivered during inhalation..

Application of a mask to a patient's face adds a certain amount of deadspace to the oro-nasal profile. Since some exhaled air remains in the area between the mask and the patient's face, applying masks with limited venting (especially without PAP pressure applied) can lead to  $CO_2$  rebreathing and/or an increase in  $ETCO_2$ . The ISO 17510-2 mask standard specifically outlines what is acceptable in terms of changes in  $ETCO_2$  due to a mask being applied. Hence it was determined to utilize the test protocol outlined in the standard to note if application of the TopOx led to an increase in  $ETCO_2$  and/or rebreathed  $CO_2$ . Results showed that  $ETCO_2$  did increase when the mask was applied, largely due to the deadspace added to the oro-nasal area. However, this increase was well within the acceptable range described in the single-fault condition tests of the standard, where positive pressures are absent from the mask, as is the case when the TopOx is applied.

Application of the TopOx mask significantly reduced the liter flow requirements to achieve specific  $FiO_{2s}$  in the tested conditions, suggesting significant oxygen conservation benefits when using the TopOx mask versus standard CFO.


 <b>Valley Inspired Products</b>		<b>Title:</b> Topox Oxygen Mask Evaluation	
<b>Origination Date:</b> June 5, 2013	<b>Revision Date:</b> July 3, 2013	<b>Document Number:</b> 13002	<b>Rev:</b> A

Exceeding the volume capacity of the O<sub>2</sub> storage bag on inspiration did not significantly affect the breathing pattern at 15 or 30 BPM. There was some resistance/work of breathing when breathing in deeply with large tidal volumes (>750mL), noted by a peak-to-peak increase of about 2 cmH<sub>2</sub>O in the pressure profile, and visual inspection during inhalation identified the mask being slightly drawn in to the mannequin head, but the looseness of the one way flappers and mask fit essentially allowed enough flow to pass through the mask to the lung without causing the flow pattern to be significantly altered or abruptly halted. If the mask were to be held forcefully to the face and/or the one-way valves blocked, then the possibility of additional resistance/work of breathing exists.

It is understood that the TopOx mask tested for this evaluation was in prototype stage, so any modifications to the product going forward may impact device performance and the results presented here should not be considered applicable to any device revisions.

All testing, data analysis and initial report generation completed from May 23—June 5, 2013. Evaluation of liter flow equivalency and RevA report completed July 2-3, 2013.

Signed:




---

Ryan Diesem, Research Manager